6610069601

Venus® Diamond Flow

Venus® Diamond Flow

- DE
- Gebruiksaanwijzing
- GB
- Instructions for use
- FR
- Mode d'emploi
- ES
- Instrucciones de uso
- IT
- Istruzioni per l'uso
- PT
- Instruções de uso

- NL
- Gebruiksaanwijzing
- SE
- Bruksanvisning
- DK
- Brugervejledning
- NO
- Bruksanvisning
- FI
- Käyttöohjeet

 **Manufacturer:**
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
Made in Germany

Distributed in USA/ Canada exclusively by:
Kulzer, LLC
4315 South Lafayette Blvd.
South Bend, IN 46614-2517
1-800-431-1785

Caution: Federal law restricts this device to sale by or on the order of a dentist.

Dentistry – Polymer-based restorative materials

ISO 4049 (type 1, class 2, group 1)

Explanation of symbols on labelling	
	Medical device
	Manufacturer
	Batch code
	Use-by date
	Date of manufacture
	Consult instructions for use
	Storage temperature
	Keep dry
	Keep away from sunlight / radioactive sources
	Contains hazardous substances
	Contains nano materials
	Unique Device Identification

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Instructions for use

GB

Medical Device – for use by dental health care professionals only.

Intended purpose

Venus Diamond Flow is a flowable, light-curing, radiopaque nano-hybrid composite. It is indicated for adhesively bonded, tooth colored restorations of anterior and posterior teeth including occlusal surfaces. Its radiopacity is equivalent to 3.0–3.5 mm of aluminium (dentine = 1 mm, enamel = 2 mm).

Indications

- Enlarged fissure sealing
- Cavity lining – as the first layer for Class I and II cavities
- Class V fillings
- Minimally invasive Class I and II fillings in areas not subjected to masticatory forces
- Minimally invasive Class III fillings
- Small repairs of direct and indirect restorations combined with a suitable bonding agent
- Splinting of mobile teeth

Contraindications

Application to dentine close to the pulp (less than 1 mm). In case of known or suspected allergy to components of the product, its use is contraindicated. Do not use this product in case of known or suspected allergies against (meth)acrylate compounds.

Patient target group

This dental material is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user

For use by dental health care professionals only.

Clinical benefit

Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthodontically.

Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases.

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Composition

Barium Aluminium Boro Fluor Silicate Glass, E4BADMA, UDMA, Ytterbiumfluorid, Silica, Titanium Dioxid, BHT, Aminoben-zoacridester, Camphorquione, Fluorescent Pigments, Metallic Oxide Pigments, Organic Pigments, Venus Diamond Flow contains approximately 41% filler by volume, its filler particle size is 20 µm – 5 µm.

Shades

2) Pure Shades for monochromal layering
Corresponding VITA® shades

LIGHT	A1, B1, B2, C1
MEDIUM	A2, A3, B3, B4, C2, D2, D3
DARK	A3.5, A4, C3, C4, D4
BLEACH	BXL, BL

2) AT Shades for monochromal and polychromatic layering

Opaque	Baseliner, OM (Opaque Medium)
Universal	BXL (Bleach Extra Light), A1, A2, A3, A3.5, A4, B2, HK A2.5, HK A5
Translucent	CL (Clear)

Kind of application

Shade selection should be done using the Venus Diamond shade guide before inserting the dental dam as the teeth change their colour noticeably once dry.

1. For ideal results, use Ivory rubber dam (or similar) to control contamination from moisture, blood and saliva.
2. Remove any caries and protect any hard dental structures according to the general guidelines for adhesive fillings.
3. Then clean and cure the cavity.
4. A quality adhesive system (e.g., BOND® Bondings by Kulzer) should be applied to the enamel and dentine according to the manufacturer's instructions before inserting the filling material. Adhesive-treated cavities must not be contam-inated with moisture, blood or saliva.
4. Apply Venus Diamond Flow in thin layers (max. 2 mm, baseliner max. 1 mm) to the cavity. When applying the Venus Diamond Flow, ensure that it flows slowly and evenly to prevent the formation of air bubbles. If necessary, remove any air bubbles using a probe before polymerisation.
5. Polymerise each layer according to the following table using an suitable blue light curing unit (wavelength peak at 450–480 nm; light output of 1550–600 mW/cm²; e.g. Kulzer Translux®). If the tooth surface being treated is not to be instrumented by a diamond stone, carbide bur, drill or other abrasive instrument, the surface should be cleaned with pumice or a fluoride free polishing paste before restoration is placed. A non-sticking placement instrument is recommended for use in clinically placing, sculpting and shaping the material.

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Recommended Venus Diamond Flow polymerisation times using Translux curing light or similar halogen curing lights:

Shade	Time	Curing depth
A1, A2, A3, A3.5, A4, B2, OM, BXL, CL, HK A2.5, HK A5, LIGHT, MEDIUM, DARK, BLEACH	20s	2 mm
Baseliner	20s	1 mm

6. During polymerisation a dispersion layer forms on the surface. This layer should not be touched or removed if further composite layers are to be applied.

7. Venus Diamond Flow can be prepared and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishes (e.g. Venus Supra) and polishing brushes are suitable for preparation.

Instructions

During long restorative procedures the surgery light should be temporarily directed away from the treatment area to prevent premature initiation of the polymerisation reaction. The polymerisation times specified above are minimum times and should not be reduced so that an optimal polymerisation rate, and the associated physical material properties, can be reached. Afterward, apply and cure the flow composite. After ejection of the material from the syringe, release the plunger to prevent continued flow. At the end of treatment, remove the hygienic protective cover from the syringe and discard it. Discard the used application cannula and carefully reseat the syringe with the cap provided. Do not reuse the syringes if there is evidence of contamination with body fluids.

Storage conditions

Store at temperatures of -10 to 25°C (14 to 77°F). Do not expose to direct sunlight. Protect from ionising radiation. Store in a dry place. Do not use the material after the expiration date.

Warnings and safety instructions

May cause skin and eye irritation. Avoid contact with skin. Skin contact may cause sensitization. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist. May damage fertility or the unborn child.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

Application notes

Visually inspect for damage before use. Damaged products must not be used. Composites should only be used at room temperature. In the case of deep cavity preparations in close proximity to the pulp or a direct pulpal exposure, cover the

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Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG