3M ESPE

RelyX™ Unicum 2 Automix

Self-adhesive Resin Cement

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ENGLISH

Product Description
RelyX™ Unicum 2 Automix is a dual-curing, self-adhesive resin cement supplied in an automix syringe. It is used for the adhesive cementation of indirect all-ceramic, composite, or metal restorations and for posts and screws. When RelyX Unicum 2 Automix is used, bonding and conditioning of the tooth structure are not necessary.

Unlike other cementation materials, RelyX Unicum 2 Automix is characterized by high stability in combination with good flowability under pressure (structural viscosity). The cement is available in various shades.

RelyX Unicum 2 Automix contains bis-functional (meth)acrylates. The proportion of opaque fillers is about 43% by volume; the grain size (D90%) is about 12.5 µm. The mixing ratio, based on volume, is 1:1 part base paste : 1 part catalyst.

These Instructions for Use must be kept for reference for the duration of product use. Please see the pertinent information for use for details on all of the products mentioned below.

Indications

- Final cementing of all-ceramic, composite, or metal inlays, onlays, crowns and bridges; 2-3-unit inlay/overlay bridges and 3-unit inlay/onlay bridges (excluded for patients with bruxism or paraodontitis)
- Final cementing of all-ceramic, composite, or metal restorations on implant abutments

Indications such as Maryland and inlay/onlay bridges (resin bonded bridges) demand especially high adhesive bond strength. Regardless of the manufacturer of the cement and restoration, these indications may be exposed to a higher risk of decementation. To achieve an optimal result when using RelyX Unicum 2 Automix, please refer to the sections "Tooth Preparation" and "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Precautionary Measures

For Patients and Dental Personnel

Base paste: contact with eyes may cause severe eye damage. Wear eye protection to prevent injury. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Base paste: contains irritants which may trigger an allergic respiratory reaction in certain individuals. This product may not be appropriate for use in those individuals with known sensitivity to sulfites since a cross-reaction may occur with sodium persulfate.

For Patients

- This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate and/or peroxide allergies.
- If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use.

For Dental Personnel

- This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash with water.
- The use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If the product contacts the glove, remove and discard the glove, wash the hands immediately with soap and water and then re-glove.
- If an allergic reaction occurs, seek medical attention as needed.

MSDS can be obtained from www.mmm.com, or contact your local subsidiary.

Precautions During Processing

- Seal temporary restorations using a eugenol-free product (e.g., with RelyX™ Temp NE). The use of temporary cements containing eugenol can inhibit the polymerization process of RelyX Unicum 2 Automix during the final cementation. Do not use hydrogen peroxide (H₂O₂) as the residues can adversely affect the adhesive strength and setting reaction of RelyX Unicum 2 Automix.

Pulp Protection

To avoid pulp irritation, cover areas in close proximity to the pulp by applying small amounts of hard-setting, zinc phosphate material. The use of zinc phosphate during the cementing of the later restoration, the pulp protection should be applied prior to taking an impression for the final restoration.

Tooth Preparation

Preparation of the Cavity/Tooth Stump

- Prior to final cementation, clean the prepared stump or the cavity thoroughly with pumice slurry, rinse with a water spray, and lightly dry in only 2-3 bursts of water-free and oil-free air, or use cotton pellets to try it off. Do not overdry!
- The cavity should be just dry enough that the surface has a slightly glossy appearance. As is the case with any permanent cement, over drying can lead to post-operative sensitivity.

Do not use substances such as desensitizers, disinfectants, astrigenes, dentin sealants, rinsing solutions containing EDTA, etc., at the final cleaning with pumice slurry and water. Their residues may have a detrimental effect on the bonding strength and setting reaction of the cement.

Preparation of Root Canals

Use the root canal as airtightness as usual (root canal filling with gutta-percha and removal of the root canal filling, leaving 4 mm of gutta-percha apically). Clean the root canal with a 2.5%-5.25% solution of sodium hypochlorite (NaOCl). Rinse immediately with water and dry with paper points.

- Discard immediately with water and dry with paper points.

We recommend the use of a rubber dam during the cementation of posts.

Preparation for Maryland and Inlay/Onlay Bridges

Abutment teeth must have an adequate enamel surface for bonding. They should be healthy or only slightly restored and the periodontal conditions should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- Prepare retentive elements such as circular rests and/or approximal grooves.
- Use a rubber dam and keep the cavity free of any contamination during cementation.

Etch the enamel surface of the cavity with 37 % phosphoric acid for 15–20 sec. Then rinse thoroughly with water and dry with water-free and oil-free air. In the case of exposed dentin, make sure to avoid the inclusion of the enamel to avoid post-operative sensitivity.

Preparatory Measures

- Remove the temporary restoration and thoroughly clean the tooth surface of any residue from the temporary cement.

- Make a trial insertion of the final restoration and check its fit and contact points.
- If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely.

For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.

- Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment and until the final cementation.

Pre-treatment of Restorations

Pre-treatment of Metal Restorations

Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- Blast the restoration surface to be luted with aluminum oxide ≤40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Etchable Glass Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.
- Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use, e.g.: – RelyX™ Ceramic Primer; allow to react for 5 sec, then air dry so that the solvent evaporates completely.

Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Blast the restoration surface to be luted with aluminum oxide ≤40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Alternative 1:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet™ Prep and the blast-coating agent CoJet™ Sand from a distance of 2–10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand, respectively.
- Blow away any residues of the blasting agent with water-free and oil-free air.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Composite Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Blast the composite restoration surface to be luted with aluminum oxide ≤40 µm or arfinium.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Pre-treatment of Glass Fiber-reinforced Posts

Please follow the instructions for use of the post to be used. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Clean the post with alcohol and dry it with water-free and oil-free air.
- When using RelyX™ Fiber Post application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with the relevant instructions for use.

Pre-treatment of Implant Abutments

Please follow the recommendations for the relevant restorative material (e.g., metal, zirconia).

Pre-treatment of Maryland and Inlay/Onlay Bridges

Metal, Zirconia and Aluminum Oxide Surfaces:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet™ Prep and the blast-coating agent CoJet™ Sand from a distance of 2–10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For etchable glass ceramics, please refer to "Pre-treatment of Etchable Glass Ceramic Restorations".

Times

The processing and setting times depend on the ambient and oral temperature. The times shown are based on conditions relevant for practice. As is the case with every composite cement, the setting of RelyX Unicum 2 Automix cement slows down significantly at room temperature.

RelyX Unicum 2 Automix is a dual-cure cement and therefore also sensitive to natural or artificial light. The working time is significantly reduced during application under operating lights!

	min:sec
Working time from start of mixing:	02:00
Light curing time:	
– single surface, from occlusal	00:20
– any other surface additional	00:20
RelyX Fiber Post posts, from occlusal	00:40
Self-curing:	
Polymerization start after start of mixing	02:30
Polymerization end after start of mixing	06:00

Application

- Select a RelyX Unicum 2 Automix syringe containing the desired color and take it ready for use.
- If allergic reactions occur, seek medical attention as needed.
- Remove the automix syringe from the foil bag and discard the bag.
- Notes down the date the syringe was removed on the syringe label.
- Remove and discard the sealing cap from the automix syringe. Check the syringe openings for blockage; remove any paste plugs.
- **The sealing cap must not be reused to close the syringe.** Because doing so might introduce air bubbles into the material.
- Squeeze out a small quantity of paste onto a block to equal the base paste and catalyst paste in the automix syringe. Discard the paste which has been squeezed out.

During every application:

Protect the working area from contamination with water, blood, saliva and sulcal fluid during the application and setting phase.

- Remove and discard the mixing tip remaining on the automix syringe from the previous application.
- Attach one of the two types of mixing tips ("Standard" or "Wide") and secure it by turning it to the right.

- The mixing tip "Standard" is used without an Intraoral Tip for application in easily accessible areas.
- The mixing tip "Wide" is used for the application with the Intraoral Tip or Endo Tip.
- When using the mixing tip "Wide," attach an Intraoral Tip or an Endo Tip.
- Squeeze out and discard a peppercorn-size quantity of RelyX Unicum 2 Automix until an evenly mixed paste in a homogeneous color flows out of the tip.

The paste requires a certain amount of time to flow through the mixing tip. The flowthrough speed cannot be accelerated by increasing the pressure on the plunger.

As soon as the pressure on the plunger decreases, the material flow stops and the paste begins to set up. Do not use force to press out paste stops and the paste begins to set up. Do not use force to press out paste stops and the paste begins to set up. Do not use force to press out paste stops and the paste begins to set up.

- Apply RelyX Unicum 2 Automix evenly to the entire cavity and as appropriate to the bottom side of the inlay/onlay, or if the crown with cement.
- Keep the opening of the mixing tip "Standard" or of the Intraoral or Endo Tip immersed in the restoration during the entire application to prevent the inclusion of any air bubbles.

- Seat the restoration firmly and stabilize long enough for the cement to set fully.
- Leave the used mixing tip on the automix syringe as a cap until the next application.

Application in the Root Canal

See instructions under "Application" for handling of the automix syringe. Do not use Lentulo-Spirale to insert the cement in the root canal as this can excessively accelerate setting.

- Attach an Endo Tip to the mixing tip "Wide" for application in the root canal.
- Insert the Endo Tip as deeply as possible in the root canal and apply RelyX Unicum 2 Automix, beginning apically. Keep the tip of the Endo Tip immersed in the cement and slowly move the Endo Tip upwards as the level of the paste rises.

- Do not remove the Endo Tip from the cement until the root canal has been completely filled.
- This so-called irrigation filling should not be carried out in less than 5 sec; this will minimize the entrapment of air bubbles.

- Place the post in the root canal filled with cement; apply moderate pressure to hold it in position. We recommend rotating the post slightly during insertion to avoid the inclusion of air bubbles.

Removal of Excess

Notes on oxygen inhibition: As is the case with all composite materials, RelyX Unicum 2 Automix is subject to oxygen inhibition, i.e., the upper layer (about 50 µm) which is in contact with atmospheric oxygen during the polymerization does not harden. Less sufficient curing permits the removal of the uncured layer during shaping/polishing without leaving behind a deficit.

From Restorations

- Using light-curing:
 - After brief curing (about 2 sec using a standard polymerization device), use a suitable instrument (e.g., scaler) to remove the excess. Larger volume excess can be removed more easily.
 - Immediately light-cure after removal (see Polymerization and Shaping).

Using self-curing:

- Use a suitable instrument (e.g., scaler) to remove the excess during the self-curing phase (about 3 min after beginning the mixing in the "gel phase").
- Immediately cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.

Alternative 1:

- Remove the excess before polymerization, e.g., using a sponge pellet, while using a suitable instrument to hold the restoration in position.
- Cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.
- Immediately light-cure or wait for the self-curing to finish (see Polymerization and Shaping).

From Posts

- Remove the cement with a suitable instrument or a cotton pellet.

Polymerization and Shaping

- We recommend light-curing the cement through the restoration when doing ceramic and composite work. Select the exposure times appropriate for the number of surfaces (please refer to "Times"). Polymerize the cement through the post when using translucent posts. The exposure time depends on the transparency of the post being used; for RelyX Fiber Post, it is 40 sec.
- Shape any remaining uneven areas on the edges of restorations and polish the marginal area with diamond devices, aluminum oxide coated discs (e.g., SoF-Lex™), and diamond polishing paste.
- Then check the occlusion.

After Shaping

- Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing on implant abutments.

Hygiene and Disinfection

- We recommend the use of commonly available hygienic protective covering to avoid any contamination of the RelyX Unicum 2 Automix syringes during the treatment.
- Clean contaminated RelyX Unicum 2 Automix syringes with cleaning agents commonly used in the dental practice.
- Use a cloth soaked in a disinfectant commonly used in the dental practice to disinfect the automix syringe.

Notes

- RelyX Unicum 2 Automix can prematurely polymerize when exposed to natural or artificial lighting. Therefore avoid intensive light exposure during application.

Storage and Stability

Store RelyX Unicum 2 Automix in the pouch at 15°–25 °C/59°–77 °F. After removal from the pouch, use the product within 6 months and before the lapse of the expiration date. Constantly high humidity accelerates the setting and must be avoided.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for the intended application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty contract, negligence, or strict liability.

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3M, ESPE, CoJet, Lava, RelyX and SoF-Lex are trademarks of 3M Company or 3M Deutschland GmbH.
Used under license in Canada. © 2012, 3M. All rights reserved.