



perio maintenance BDR



Biofilm disruption



TK1-1S
F01001



Short probe

Graduated every 3mm, for examining shallow and medium pockets (< 4mm) and for the maintenance of simple cases.



TK1-1L
F01004



Long probe

Examination and maintenance of medium to deep pockets (> 4mm).
Diagnosis aid during the debridement and irrigation of pockets.

The TK1 probe tips are used without pressure following the contour of the pockets and skimming over the root surface.



TK2-1L
F02162



Maintenance of the premolar and molar sectors, left-oriented

Maintenance of moderate to deep pockets and furcations.
Equivalent to the Nabers probe.



TK2-1R
F02161



Maintenance of the premolar and molar sectors, right-oriented

Complementary to the TK2-1L tip for the maintenance of moderate to deep pockets and furcations.
Equivalent to the Nabers probe.