



# Instruments for Injection



## PAROJECT®


Injection syringe for intraligamental local anaesthesia.

PAROJECT® is an unthreatening pen style intraligamental syringe with a controlled dosage permitting a slow non-traumatic injection.

- Compact and user-friendly (delicate pen-type design)
- Unthreatening in appearance
- Injects only 0.06 ml. per "click", permitting a slow non-traumatic injection
- Made from high-grade stainless steel
- Can be sterilized by any method
- Material and construction guaranteed for 5 years

### PerioDontal Ligament Anaesthesia PDLA

After a short while of practising the intraligamental technique in general facilitates everyday treatments, and saves time. The patient will experience that dental local anaesthesia no longer has to be tantamount to a numb cheek, tongue or lip for several hours. The intraligamental technique is a gain to the dentist as well as to the patient.



### Advantages and Indications

- Immediate onset.
- Profound single-tooth anaesthesia in connection with e.g. simple extractions, cavity preparations, endodontic treatments, etc.
- Overcomes failed conventional techniques - an efficient supplement to insufficient blocks and infiltrations.
- Teeth in more than one quadrant can be treated at the same sitting.
- Gentle pain control due to a very small injection dose.
- No soft tissue numbness.
- Replacement to most palatal injections.
- Substantially reduces the risk of haematoma - important when treating haemophilic patients (where conventional injection techniques might be contra-indicated).



**RØNVIG**  
Dental Mfg. A/S

Gl. Vejlevej 59 · DK-8721 Dagaard · Denmark  
Tel.: (+45) 70 23 34 11 · Fax: (+45) 76 74 07 98  
e-mail: [export@ronvig.com](mailto:export@ronvig.com) · [www.ronvig.com](http://www.ronvig.com)