

# MicroCab+™


## Self-Contained Dust Cabinet

**MicroCab+** is the most effective sandblast and grinding containment system available. Offering a large internal work space, the MicroCab+ still maintains a small footprint. Optional wall mount bracket makes the unit even more versatile without taking up valuable counter space. Custom ergonomic design offers integrated handle, bright fluorescent light and ultra quiet fan/filter. The MicroCab+ is a must have in any lab or operatory. 110 volt (220v optional).



**MicroCab+**  
Includes 5 windows  
6.2 (h) x 8.7 (w) x 10.5 (d)  
Ref # 91010 (120v)  
Ref # 91354 (220v)

**Disposable windows**  
Ref # 91286 10/pk  
**Wall mount bracket**  
Ref # 91350



**MacroCab**  
7 (h) x 10 (w) x 14 (d)  
Ref # 58600

# MacroCab™

## Self-Contained Dust Cabinet

Compact, yet larger, unit offering an efficient and comfortable design where you only insert your fingers with handpiece to work. Bright, fluorescent light provides ideal illumination. Featuring built-in dust collection/filtration system. It is easy to clean and the filter lasts for years. It is quiet, made of heavy gauge aluminum with a baked enamel finish. 110 volt (220v optional).



**Safety Shield (10 pack)**  
MacroCab Ref # 201360  
MicroCab Ref # 201359

# Safety Shields

## Disposable Window Shields

Extend the life of your MacroCab or MicroCab+ window with disposable window shields. When placed under cabinets original window, these safety shields protect it from harsh abrasive.

[www.danvillematerials.com](http://www.danvillematerials.com)

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance.

92195 Rev A