

endodontics



Retreatment and obturation



Long retreatment tip, 40mm, taper 4%

Rapid removal of broken instruments in the middle third of wide, straight canals.



Long retreatment tip, 40mm, diamond-coated 30 μ m, taper 4%


Retreatment of very hard material in the middle third.



Long Titanium-Niobium tip, 25mm, taper 3%

Retreatment in the apical third and long, straight canals.

ET25 tips can be pre-formed for the treatment of curved canals.



Fine condenser, length 40mm, taper 4%

Lateral condensation of gutta percha by heating effect. It is used dry, without irrigation.



ET40
F88012



ET40D
F88014



ET25L
F88022



SO4
F88009