



CLEANING PERFORMANCE: “EC.WASHTEST” CLEANING PERFORMANCE TEST/WASHING INDICATOR

- Chemical indicator evaluates the effectiveness of thermo-disinfectors and controls the cleaning conditions in the automatic washing machines.
- The non-toxic chemical reagent dissolves with most of the cleaning agent and is swept along with the washing fluid.
- 250 pcs/box.

CLEANING PERFORMANCE: “WT HOLDER” WASH TEST HOLDER

- Holder for Wash Test
- Made of stainless steel
- Highly resistant to corrosion and abrasion


CLEANING PERFORMANCE: “ULTRASONIC CHECK” ULTRASONIC CLEANING INDICATOR

- Ultrasonic cleaning processes (cavitation) allow disinfection even in areas difficult to access by other cleaning methods. This indicator verifies the effectiveness of the ultrasonic washing equipment.
- Ready to use either with or without detergents and cleaning solutions.
- Insert in the machine for at least 3 minutes between 18-70 ° C.
- Extremely reliable and stable product.
- Rapid read-out.
- 30 pcs/box.

CLEANING PERFORMANCE: “PROTEIN CONTROL” CONTROL FOR SURFACE PROTEIN

- An accurate disinfection ensures a reduction in bacterial growth and resulting nosocomial infections.
- This system is able to detect any proteins and allergens on the surfaces after the cleaning process.
- High-reliable and high-sensitive product.
- Intuitive interpretation of the result.

REF. EC.WASHTEST



REF. ULTRASONIC CHECK



REF. PROTEIN CONTROL

