

AccuFilm® FASTCHECK

Occlusal Articulating Film

Articulation. Simplified.

- ▶ Works for a single-tooth or entire quadrant.
- ▶ Built-in holders eliminate the need for forceps.
- ▶ Consistently produces sharp, easy-to-interpret marks.
- ▶ Easily identify high contacts making for easier adjustments.
- ▶ Simultaneous marking of both arches with its double-sided film.
- ▶ Incredibly thin film for accurate exposure of high occlusal contacts.
- ▶ For use on wet or dry surfaces.



AccuFilm®, the world-renowned brand recognized for accurately marking occlusion for decades introduces its simplest film to date! AccuFilm® FastCheck incorporates the accuracy that dental professionals have come to rely on with the convenience of built-in handles.



Film strip easily detaches from protective backing by using the convenient handles and gently pulling in opposite directions.



19 MICRONS ✓
FILM THICKNESS

SKU # PRODUCT DESCRIPTION

S054	AccuFilm® FastCheck Double-Sided Articulating Film — Blue/Blue (100) Strips
------	---

See next page for our complete line of AccuFilm Articulating products.

To learn more about AccuFilm FastCheck and our complete line of AccuFilm articulating products, contact your Parkell Representative or go to parkell.com.

AccuFilm® Booklets

SKU # PRODUCT DESCRIPTION

S050	 AccuFilm® Arch Double-Sided Marking Paper Booklet—Red/Blue (6) booklets (12 sheets/booklet, 72 sheets ttl.).	
S053	 AccuFilm® Double-Sided Marking Film Booklet—Red/Black (5) booklets (50 sheets/booklet, 250 sheets ttl.).	

Double-sided film thickness is 19 microns.

AccuFilm® Pre-Cut Film Canisters

SKU # PRODUCT DESCRIPTION

S017	 AccuFilm® Double-Sided Marking Film—Red/Black (1) canister of 280 pre-cut strips.	
S018	 AccuFilm® Double-Sided Marking Film—Black/Black (1) canister of 280 pre-cut strips.	
S021	 AccuFilm® Double-Sided Marking Film—Red/Red (1) canister of 280 pre-cut strips.	
S014	 AccuFilm® Single-Sided Marking Film—Black (1) canister of 280 pre-cut strips.	
S015	 AccuFilm® Single-Sided Marking Film—Red (1) canister of 280 pre-cut strips.	

Single and double-sided film thickness is 19 microns.

AccuFilm® IV Brush-On Liquid

SKU # PRODUCT DESCRIPTION

S020	 AccuFilm® IV Brush-On Marking Liquid (1) bottle of brush-on marking liquid, (1) bottle of solvent. <i>Made in Germany.</i>	
------	---	---