

Root planing



H1

Anterior tooth root planing, diamond-coated tip 30 µm

- Diamond-coated mini-tip for simple cases in the cervical area.
- Also effective for the withdrawal of granulation tissue.

This tip should be used without pressure and above the epithelial attachment because it is abrasive.



H2L

Root planing of the premolar and molar sectors, left-oriented, diamond-coated tip 30 µm

Diamond-coated micro-probe for the treatment of furcations and narrow spaces.



H2R

Root planing of the premolar and molar sectors, right-oriented, diamond-coated tip 30 µm

Diamond-coated micro-probe for the treatment of furcations and narrow spaces.



H1
F00366



H2L
F00367



H2R
F00368