# **ENBIO**<sup>®</sup>

World's fastest B class air-removal (pre-vacuum) autoclave.







Completes of up to 50 full cycles per day.

Brings the tools ready when you need them.

Reduces the number of expensive tools needed.





Sterilizes wrapped and unwrapped tools.

Ultra compact, with fast turnaround.

Incredibly intuitive: just touch the screen.

B class air-removal (pre-vacuum) autoclave





Safety takes the stage.

Fits your working environment.

Works so silent that you will barely notice it in the room.

Shows your patients that you care.

### enblo.







1-click process

**Quality tested** 

Class B (Vacuum)



We get you covered.



### Warranty

Your peace of mind is what matters the most. We think so.



### **Backup**

You always get a replacement device for the time of service.

### **Program**

## 7 min

## **Total process time**Use for

Solid objects, non-porous objects, simple instruments (such as sissors, handles, pliers, chisels, probes, etc.), and dental handpieces; unwarpped for immediate use only.

134°C

**FAST** 

Drying

## 15 min

Solid objects, small porous objects, simple objects recessed, narrow clearance items, dental handpieces and textiles; wrapped and unwrapped.

134°C

YES

## 121°C

## 30 min

Solid objects, small porous objects, simple objects recessed, narrow clearance items, dental handpieces textiles and plastics; wrapped and unwrapped.

YES





### **Technical data**

Power supply
Installed power
Max. power consumption
Working pressure
Max. pressure
Max. process tempreature
Process chamber volume

220-230V AC/ 50/60Hz 2250 W

2250 W 15A 2,1 Bar 2,3 Bar 131°C 2,71 Max. packeges amount (circa)
Max. batch weight
Mass
Level of security
Data archiving
Water conductivity\*

\* demineralized or distilled water only

8 0,5kg 15 kg IP20 USB Me mory stick <15 µS/cm

#### Certifications

In accordance with the requirments of the EN-ISO 13060 standard and 93/42/ EEC Directive. Medical device class Ilb.